

PUETTERPRO 2 2KOMP SYS10CM

PÜTTERPRO 2 2-Komponenten-Kompressionssystem 10 cm

PZN: 11482930

Inhalt: 2 St

Hersteller: Paul Hartmann AG

Zweckbestimmung

PütterPro 2 is indicated for the treatment of chronic venous insufficiency (CVI) with and without ulcer cruris venosum (UCV) as well as chronic oedemas and an ABPI of 0,8-1,3.

Zusammensetzung

Layer 1: 92% polyester + 8% elastane, with special latex-free coating on one side; non-woven material with elastane threads; tearable; cut edges. Layer 2: Compression bandage, flesh-colored, 84% polypropylene + 16% elastane, with special latex-free coating on both sides; non-woven material with elastane threads; tearable; cut edges