

VERLAENG M V0.9X1.9MM 300

VERLÄNGERUNGSLEITUNG 0,9x1,9x300mm Schiebeklemme

PZN: 17400620

Inhalt: 1 St

Hersteller: B. Braun Melsungen AG

Zweckbestimmung

Extension Lines are used in all medical disciplines where IV therapy is performed. IV therapy is performed for :- Fluid volume maintenance,- Fluid volume replacement,- Delivery of chemotherapeutic agents, nutritional support, antibiotics and other medications. The use of the Extension Lines with IV-Solutions, which have a proven incompatibility to one of the used materials, is contraindicated

Zusammensetzung

PE-HD; PVC; PVC-P; PC; PE-LD